KaVo Cellular Optic High-Speed Handpieces

GENTLEforce LUX 3 6000B • Super-Torque LUX 3 647B • Mira LUX 3 635B





KaVo KaVo KaVo KaVo KaVo KaVo KaVo YOU enduring nature. And born to lead. We're in this together. We don't take our premier position with you for granted. Like you, we're both artist and engineer. And our latest high-speed isn't just a handpiece. It's our latest advance designed to give you another advantage. Balance. Beauty. Brute strength. A precision instrument that blurs the distinctions. A dedication to you and dentistry and driven to stay the #1 choice of dentists worldwide.



- Quietest KaVo high-speed handpiece ever
- 4-port spray for better visibility and cooling
- New hygienic grip improved tactile feedback
- 3rd generation cellular optic glass rod 5–year warranty
- Proprietary replaceable microfilter better spray quality
- Unique automatic drive air adjustment longer turbine life
- Non-return system prevents debris from entering the turbine
- Sterilizable up to 135°C
- 18-month warranty including turbine*

*If maintained with KaVo QUATTROcare (standard warranty 12 months)

Switch to KaVo — KaVo handpieces easily attach to virtually every existing tubing and light source system with a MULTIflex coupler.




You know you really want KaVo.

In study after study*, when dentists are asked what they look for most in a high-speed handpiece, they list:

- Perfect concentricity
- Powerful torque
- Durability even after repeated sterilization

KaVo high-speed handpieces bring dentists' most-requested features to the table, and even more. The KaVo reputation rests on a long history of engineering excellence and an understanding of what dentists demand. So, it's no wonder that when dentists are asked, "If you could own any handpiece, which would you choose?" More dentists choose KaVo than any other brand.

And, now there's more reason than ever to go with the leader. Our newly designed GENTLEforce LUX 6000B blows away the competition with its easy-access design, smaller head and absolutely the highest torque of *any* multiple spray port high-speed handpiece available today.



6000B Power/Torque Relations

Exceed your highest expectations

KaVo is the most desired handpiece for good reason. Each and every handpiece undergoes exhaustive quality control checks and testing before it leaves the factory. KaVo's stringent specifications and a zero tolerance for defects guarantee optimum performance with a minimum of maintenance.

Choose the newly designed 6000B or one of our standard high-speeds.

KaVo GENTLEforce 6000B - New Small head

Ultimate torque
4-port spray
For gross reduction and prep, pedo
and difficult-to-access areas

KaVo Super-torque 647B

Standard head Very high torque 3-port spray For gross reduction and prep






KaVo Mira LUX 3 635B


Miniature head Highest torque in a miniature head 1-port spray

For pedo and difficult to access areas



Clean feeling of luxury
Unique surface technology: The new
6000B features a unique Plasmatec surface
that actually works with the innate gripping
power of the glove for slip-free, comfortable
control. And since there are no grooves or
edges to harbor debris, the 6000B
provides superior hygiene.

The new angle combination between handle and head increases freedom of movement by 3° . You have to experience it to feel the difference.



New Design—smaller head, new angles and better access

Durability, consistent feel and precision balance are KaVo hallmarks. The new 6000B pushes our tradition of function and design into the 21st century.

Precise preparations are a result of your skills, concentric turbine operation, instrument angulation and head size. The ingenious KaVo GENTLEforce 6000B creates a new feeling of more space in the oral cavity when you put our new design to work.

Freedom of Movement

- · New handpiece angles for improved anterior and posterior access
- Smaller head for better vision and precision
- You work closer to the preparation
- Greater mobility even where space is limited
- Faster and easier preparations
- Gentler treatment

Exceptional Comfort and Convenience

- Quietest KaVo high-speed ever
- Improved grip for greater tactile feedback
- · Perfect balance
- Less stress and fatigue for you

Conventional high-speed

KaVo GENTLEforce LUX 6000B

All KaVo High-Speed Handpieces Feature:

Precise. Powerful. Performance.

- Optimum balance between speed and torque for faster, more consistent cutting
- Highest concentricity for smooth preparations and better defined margins
- Proprietary automatic pressure regulator ensures optimum performance and maximum turbine durability
- Durable push-button auto chuck saves time and requires minimum maintenance
- Fine spray mist for better visibility

Unmatched Flexibility

- MULTIflex Coupler for single, secure snap-on, quick connect for air, water and light
- 360° rotation due to the MULTIflex Coupler
- Absolute flexibility for diverse clinical applications with MULTIflex:
 - Air abrasion
 - Air polishing
 - Crown and Bridge removal
 - Air scaling
 - Light probe

The smaller head allows you better access to the posterior region.



New head and handle angles facilitate overall access and maneuverability.



Dentition tends to impede the freedom of movement of standard handpieces.



- Provide faster cutting 1
- Keep tooth structure cooled even in deep distal preparations
- Extend bur lifeBetter patient comfort

Outstanding Visibility
KaVo's 3rd generation cellular optic glass rods
are 25% brighter than conventional optics,
and are backed by a 5-year warranty.



¹ SIEGEL, S., Von Fraunhofer, A. The Effect of handpiece spray patterns on cutting efficiency. JADA (Vol. 133, 02/02) p.184f.

Extraordinary Durability

The results of the KaVo Quality Assurance Program are the foundation of KaVo's reputation as the handpiece market leader. Clinical evaluations have determined that KaVo handpieces offer superior performance even after repeated sterilization.

Third-Generation Optics

While others are still using first-generation lights, KaVo cellular glass rod optics represent third-generation technology that is light years ahead of anything else. In fact, KaVo was the first to introduce glass rod technology back in 1989. And since then, we've refined this technology to state-of-the-optical-art. Powerful 25,000lux is color illumination balanced for the best definition in the oral cavity. Studies show virtually no degradation of light intensity after thousands of sterilization cycles.

Extended Life

- KaVo high-speeds are proven superior after repeated sterilization cycles
- Galvanized surface coating with a durable nickel silver alloy on all models assures long life
- Plasmatec 6000B surface coating over the grip area lasts and lasts
- KaVo German precision engineering with exact tolerances to run more smoothly over its extended working life
- Individually balanced factory turbines run with absolute concentricity so you can sculpt preparations with more precise architecture, while minimizing possible microfractures

KaVo doesn't just talk durability. We back it up with an 18-month warranty on new 6000B handpieces.* The 647B and 635B feature a 12-month warranty.

Simple, easy to maintain and ultra-durable!

All KaVo handpieces are autoclavable at 135°C, and can be washed in a thermo-disinfector. Tests have shown KaVo high-speed handpieces can be sterilized thousands of times without compromising performance.



Exchangeable microfilter in spray-water channel guarantees efficient and reliable spray-cooling at all times



Anti-suck-back system due to new cartridge-design which prevents contaminated particles to enter the head-housing



Warranty terms

*The 6000B features an 18-month warranty when maintained with the KaVo QUATTROcare (standard warranty - 12-month). Genuine KaVo replacement turbines for the 6000B will also feature an 18-month (12-month) warranty.

KaVo High-Speed Handpieces	A contract to the contract to	667 E	West 1
Specifications	New KaVo GENTLEforce LUX 6000B	KaVo Super-Torque LUX 3 647B	KaVo Mira LUX 3 635B
	Ultra Torque Power	Exceptional Power	Very High Power
Specifications Specifications	#1000.3903	#0553.2090	#0553.2220
Specifications Specifications Specifications	Ideal for gross reduction and crown and bridge work, also suitable for pedo and difficult to access areas.	Ideal for gross reduction and crown and bridge work.	For pedo work and difficult to access areas
Specificatispray	4-port	3-port	1-port
Specifications Specifications Specifications	Small	Standard	Miniature
Specificonnection	Quick disconnect with MULTIflex LUX coupler		
Specifications Specifications Specifications	The highest torque ever in a KaVo high-speed. Note: head is smaller than on 649, 647 and 640 18 watts	Increased power over previous 640B 15 watts	Excellent power 13 watts
Spec Chuck System	Pc	owerful, easy to use push-button auto ch	uck
Specificat optics	Т	hird generation LUX cellular optic glass r	rod
Specifications Recommended Air Pressure	40psi/2.8 bar - minimum	36psi/2.5 bar - minimum	32psi/2.1 bar - minimum
Specifications	A unique pressure r	regulator ensures consistent cutting pow durability of the turbine	er and enhanced
Specifications Specifications Specifications Specifications Specifications Specifications Specifications Specifications	Nominated for REALITY Product of the Year for 2004 DENTAL ADVISOR ++++	Nominated for REALITY Product of the Year for 2004 REALITY Product of the Year for 2004	Nominated for REALITY Product of the Year for 2004 Nowinated for REALITY Product of the Year for 2004 Kavo Product Kavo Product ADA Product Rayland Rayland

KaVo Light sources and tubing

LS2 illuminator #0905.8500

- · Compact remote fiber optic light source
- Hanger or foot pedal installation
- · Operates up to three handpieces
- High/low intensity switch
- · Installation kit included

MULTIflex Couplers Required

1390 (only for LS2 standard tubing or other 5-hole tubing)

Tubing type to be used

· LS2 standard tubing straight or coiled use with 1390 coupler







465LRN MULTIflex LUX Coupler #0553.1550





If you already have a light source, all you need is a MULTIflex LUX coupler. Simply screw the coupler onto your existing tubing and you're ready to snap on the KaVo handpiece you've always wanted.

LCM LUX control module #0905.8400

- · Compact bulb in the coupler system
- For foot pedal and hanger installation
- Operates up to four handpieces
- · Includes module transformer, remote intensity control, installation kit and light delay system

MULTIflex Couplers Required

465LRN or 459AN (not shown)

Tubing type

LCM tubing straight 11-inch in black, gray or surf; or LCM tubing coiled 7-inch in black or gray

• 18-Month WARRANTY



When maintained with KaVo QUATTROcare, your new or KaVo factory repaired 6000B is protected with a full one-and-a-half-year warranty.

• 12-Month WARRANTY

The 647B and 635B are covered with a 1-year warranty.



· Maintenance and Service

KaVo handpieces are easy and quick to maintain with one-step KaVo Spray or the KaVo QUATTROcare. And should you need it, our legacy for support and the guaranteed factory repair service is second to none.



KaVo QUATTROcare

Clean, lubricate and purge up to four handpieces, attachments or air motors in one cycle, automatically, at the touch of a button. QUATTROcare Unit #0411.7800



KaVo Spray

One fast spray prior to sterilization is all that is needed. Rated simplest maintenance program by CRA. KaVo Spray #0411.9660



· MULITIflex nozzle for high-speed handpieces. #0411.9921



· Genuine KaVo Turbines

KaVo turbines are designed and manufactured to withstand the day-to-day repeated workload and sterilization better than any other brand. All KaVo turbines can be easily replaced in the office or sent to KaVo.


To ensure maximum performance and concentricity, bearings cannot be replaced on the 6000B.



• The best light


25,000lux, high-precision robust glass rod stands up to repeated sterilization - 5-year warranty.

· Simple, easy to maintain and ultra-durable!



WARNING!

Generic Turbines May Kill KaVo Handpieces And Threaten Your Investment



Turbines turn at more than 6,000 times per second, a tiny, imperceptible wobble or shimmy is quickly amplified resulting in poor performance and the need for replacement. That means inconvenience and higher costs.

Even if a repair service says they use genuine KaVo parts, can you really be sure? Do they electronically balance the generic turbine like KaVo does with every new cartridge? Like unbalanced tires that demonstrate faster wear and poor performance, improperly balanced turbines develop chronic wear problems that damage your handpiece, lead to faster bur breaks and prevent smooth, precise preps.



Look for the seal that certifies a genuine KaVo factory repaired handpiece. It's your ONLY guarantee that the handpiece meets original specification. It's an assurance backed by an extended one year factory warranty, the same as our new handpieces.

If you are using generic parts and non-certified repair shops, you are invalidating your warranty and jeopardizing the useful life of your handpiece investment. Most likely, you bought KaVo because it is recognized as the world's best handpiece. Why risk performance and increase costs with sub-standard parts?

Studies have revealed that dentists consider durability to be the most important handpiece feature. KaVo handpieces are the proven leader with independent evaluations that confirm superior durability and maximum handpiece life. Chances are that a problem with your KaVo handpiece is the result of faulty parts and improper repair.

A Comparison of Short and Long Term Savings			
	Generic Turbine	Genuine KaVo Turbine	
Turbine Cost	\$164¹	\$260 ²	
Expected Life	8 Months ¹	18 Months	
Average Cost Per Month	\$20.50	\$14.44	
Estimated 18-Month Cost	\$369	\$260	
Warranty Coverage	6 Months ¹	12 Months ³	

¹Average numbers received from generic turbine sales centers. Data available upon request.

²Average KaVo Turbine price (6000B – \$235, 647B and 635B – \$285 retail).

³Bigger savings with an 18-month warranty for the New GENTLEforce 6000B high-speed handpiece maintained with QUATTROcare.

Generic parts are the substantial equivalent of a disposable handpiece. If you are experiencing the following symptoms with a recently repaired handpiece, they most likely are the result of generic parts and inadequate repair.

- · Lack of bur retention leading to reduced cutting efficiency
- Non-concentric cutting resulting in less precise preparations
- · Bur breakage
- Fast wear and the need to replace parts more often
- Performance degradation after repeated sterilization This all leads to:
- · Higher repair costs
- More headaches from inconvenient breakdowns and slow service turnaround Any and all of these problems can be solved with genuine KaVo parts and certified repairs. Next time your KaVo handpiece requires maintenance and repair, do not settle for less than exceptional performance. Call KaVo today and extend the life of your investment.